Installations- und
Betriebsanleitung
Installation and
Operating Instructions
Manuel d'installation
et de maintenance

Drehschieber-Vakuumpumpen PB/ PC 0003/ 0008 A
Rotary Vane Vacuum Pumps PB/ PC 0003/ 0008 A
Pompes à Vide Rotatives à Palettes PB/ PC 0003/ 0008 A

Installations- und Betriebsanleitung PB/ PC 0003/ 0008 A
Installation and operating instructions PB/ PC 0003/ 0008 A
Manuel d'installation et de maintenance PB/ PC 0003/ 0008 A

1

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- PB 0003 A
- PC 0003 A
- PB 0008 A
- PC 0008 A

These installation and operating Instructions are valid for the following pumps:

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These vacuum pumps have been manufactured according to the latest technical standards and safety regulations. If not installed properly or not used as directed, dangerous situations or damage might occur.

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It is mandatory that these operating instructions be read and understood prior to the vacuum pump installation and start-up.

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Safety

This vacuum pumps have been manufactured according to the latest technical standards and safety regulations. If not installed properly or not used as directed, dangerous situations or damage might occur.

Application

These vacuum pumps are designed for evacuating closed vacuum chambers. The vacuum pumps can be used to evacuate air or dry gases, which are not poisonous, aggressive or explosive.

Other agents should not be transported. In case of doubt, please contact the manufacturer.

Safety advice

In this operating instruction safety measures are advised before each step. It is imperative that these safety precautions are observed.

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
Liquid and solid particles must not enter the pump.
In case of doubt consult the manufacturer.

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- 1 Gas inlet
- 2 Exhaust filter
- 3 Oil sight glass
- 4 Oil fill plug
- 5 Oil drain plug

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
- 1 Gas inlet with non-return valve
- 2 Outlet
- 3 Oil sight glass
- 4 Oil fill plug
- 5 Oil drain plug
- 6 Exhaust filter

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Transport and Packing

Minola vacuum pumps pass a rigorous operating test in the factory and are packed carefully to avoid damage in transit. The inlet flange is sealed with a plug, so no dirt can enter the pump during transport. Please check packing on delivery for transport damage.

Packing materials should be disposed of according to environmental laws or re-used.

These operating instructions are part of the consignment.

Pumps are generally shipped without oil.
Operating the vacuum pump without oil will damage the pump!

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Start-up

1. Setting-up

The pump must be set up or mounted horizontally on a flat surface.

Ambient temperature: 12 to 30° C

Ambient pressure = Atmosphere

In order to avoid overheating of the pump, an undisturbed fresh air-flow to the pump is necessary.

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2. Inlet connection

The inlet flange (Fig. 3.1.1/4.1.1) can be connected to a vacuum-tight flexible hose or pipe.

Restriction of the pipes must be avoided in order not to decrease the displacement of the pump.

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3. Oil filling

The vacuum pump always leaves the factory without oil. Operation without oil destroys the vacuum pump! Therefore oil has to be charged before start-up.

For that, open the oil fill plug (fig. 3.1.4/4.1.4) and charge the oil until the oil level reaches the upper third of the oil sight glass (fig. 3.1.3/4.1.3). After that replace the oil fill plug again.

Type and quantity of oil please see "Maintenance", pages 9 and 10.

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Conditions for installation

The pump should be installed in a position that ambient air can flow around the pump for air cooling purposes. Exhaust gas must be allowed to escape from the vicinity of the pump to avoid overheating. When installed in restricted areas a good air flow across the pump is required.

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Please note: When vacuum pump RB, RC 0003 A is switched off while under vacuum gas and oil may back flow into inlet pipe!

To avoid this a non-return valve or shut-off valve has to be installed.

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Electrical installation may only be conducted by a specialist. Regulations following EMV Directive 89/ 336/ EEC, low Voltage Directive 73/ 23/ EEC, and the appropriate EN Standards have to be applied as well as VDE/ EVU regulations and local or national regulations. The operator of the vacuum pump has to inform the manufacturer, if electric or electromagnetic interference from his mains is to be expected.

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Built-in vacuum pump:

Connect the cable with the three skinned strands according to the switch system (Fig. 7.1)

Laboratory vacuum pump:

This vacuum pump is wired ready for installation

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
and provided with a standard plug and a start-stop switch.

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Both versions are protected against thermal overloads of the motor by a thermal-lag switch. When the motor has cooled down the vacuum pump starts again automatically.

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Maintenance

The vacuum pump must be switched off and secured against accidental switch-on for all maintenance.

1. The oil level should be checked at least once a day.

2. If the oil level in the oil sight glass (fig. 3.1.3/ 4.1.3) is in the lower third, more oil must be added.

3. Changing of oil

Oil must be changed after the first 100 hours of operation. Further oil changes depend on operating conditions. The oil must be changed after 500-2000 hours of operation, but at least semi-annually. If there is considerable pollution or water in the oil sump (either as pure water or as an emulsion) oil changing is required as soon as possible.

To change the oil, the warm pump must be switched off and ventilated to reach atmospheric pressure. Drain the oil through the oil drain plug (fig. 3.1.5/4.1.5). Refasten the oil drain plug and fill fresh oil through the oil fill plug (fig. 3.1.4/ 4.1.4).

Used oil should be disposed of according to environmental laws.

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5. The exhaust filter (Fig.10.1/10.2) must be changed every 500 hours of operation. If there is considerable pollution causing increased power consumption, increased temperature and/ or decreased separation it must be changed earlier. For that unscrew the filter element (Fig. 10.1 by hand/ Fig. 10.2 by screw driver) out of the housing. After assembling fasten the new filter element hand tight.

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Amsterdam, Barcelona, Basel, Birmingham, Brussels, Copenhagen, Dublin, Gothenburg, Helsinki, Istanbul, Kuala Lumpur, Melbourne, Milan, Montreal, Moscow, New York, New Plymouth, Oslo, Paris, San Jose, Seoul, Singapore, Taipei, Tokyo, Vienna

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Dr.- Ing. K. Busch GmbH